



Connect your 3Shape Trios to Dental Art Solution Workflow

3Shape TRIOS

Send your 3Shape TRIOS scans to Dental Art Solution using 3Shape Communicate™. TRIOS order forms can be customized to match specific information requirements and offerings from our lab, making it easy to create new orders, attach 2-D images and convey all case aspects to our lab. Minutes after your scanning is complete, our lab technicians can receive your TRIOS digital impressions and immediately start designing your traditional or implant restoration. When your design is complete, a 3-D model will automatically be sent to the 3Shape Communicate site for you to review before we begin fabricating your restoration.

- Go to us.3shapecommunicate.com in a web browser .
- Connect with Dental Art Solution as Lab by searching dentalartsolution@gmail.com
- After connecting to us as a Lab, select Dental Art Solution when sending files.